



**SOPs for the operation of the capillary biofilm reactor  
available for the PROMICON consortium  
Deliverable 4.3**

31. May 2024

Hannah Berreth, Katja Bühler  
Helmholtz Centre for Environmental Research – UFZ, Leipzig, Germany

**PROMICON**

**Harnessing the power of nature through PROductive Microbial CONSortia in  
biotechnology - measure, model, master**



**Prepared under contract from the European Commission**

Grant agreement No. 101000733

EU Horizon 2020 Research and Innovation action

Project acronym: **PROMICON**  
 Project full title: **Harnessing the power of nature through PROductive Microbial CONSortia in biotechnology - measure, model, master**  
 Start of the project: June 2021  
 Duration: 48 months  
 Project coordinator: Dr. rer. nat. Jens O. Krömer  
 Helmholtz Centre for Environmental Research - UFZ  
<https://www.ufz.de/>

Deliverable title: SOPs for the operation of the capillary biofilm reactor available for the PROMICON consortium

Deliverable n°: D 4.3  
 Nature of the deliverable: Report  
 Dissemination level: Public

WP responsible: WP4  
 Lead beneficiary: UFZ

Citation: Deliverable D 4.3 EU Horizon 2020 PROMICON Project, Grant agreement No 101000733.

Due date of deliverable: Month 36  
 Actual submission date: Month 36

Deliverable status:

Version	Status	Date	Author(s)
1.0	Submitted	31.05.2024	UFZ

## Table of contents

1. Introduction.....	5
2. Methodology.....	6
2.1 Strains used in this study .....	6
2.2 Growth media used.....	6
2.3 Brief introduction to the CBR: Set-up and operation.....	9
2.4 Hydrogen analysis.....	10
2.5 CDW analysis.....	10
3. Results and discussion.....	11
3.1 Triple species biofilms comprising <i>R. palustris</i> , <i>Synechocystis</i> , and <i>P. taiwanensis</i> .....	11
3.2 Impact of nitrate on triple species biofilms including <i>Synechocystis</i> .....	12
3.3 Establishing triple species biofilms including <i>Tolypothrix</i> sp. PCC 7712 for H <sub>2</sub> production.....	13
4. Concluding remarks .....	14

## 1. Introduction

In industrial biotechnology utilizing whole-cell biocatalysts, most processes rely on mono-septic suspended cultures of the respective production strains. This leads to challenges regarding contaminations and also ongoing evolution or the formation of subpopulations, which negatively impact space-time-yields. Mutually dependent co-cultures can stabilize each other and lead to overall better production. Also, in contrast to suspended cultures it is possible to grow the organisms in a surface attached format, termed biofilm, which represent a very common microbial life style.<sup>1</sup> In biofilms organisms grow in tightly packed aggregates that are often attached to surfaces. Upon attachment these microbes start to excrete extracellular polymeric substances (EPS), mainly composed of carbohydrates, proteins and extracellular DNA, which serve as protection shield and concrete, giving biofilms a three dimensional structure.<sup>2</sup> Hence, biofilms are naturally immobilized high cell density biocatalysts that can be operated in continuous processes in which the substrate is supplied constantly to the biocatalyst.<sup>3</sup>

Here, we strive to establish an artificial microbial consortium, which grows as biofilm and is capable of producing hydrogen in a continuous mode. The most important strain in this consortium is *Rhodopseudomonas palustris* (*R. palustris*), which under redox limiting conditions (anaerobic cultivation) produces hydrogen, while growing on organic acids like acetate or lactate and thus serves as labor strain. *R. palustris* belongs to the Gram-negative purple non-sulfur bacteria and is well known for its exceptional metabolic versatility and robustness, and its ability to flexibly switch between these modes of metabolism depending on the environmental conditions.<sup>4</sup>

*R. palustris* is cultivated photoheterotrophically under anaerobic and nitrogen limiting conditions. Hence, the strain utilizes light for energy, and obtains carbon and electrons from the organic acids, while fixing nitrogen in a nitrogenase catalyzed reaction. As a side product, hydrogen is produced in reasonable amounts, depending on the specific reaction conditions.

This organism was teamed up with *Pseudomonas taiwanensis* VLB120 as a balancer strain (see also Del 3.7). *P. taiwanensis* is a chemoheterotrophic organism, which is growing aerobically and is well known for its excellent biofilm forming abilities. We hypothesized, that *P. taiwanensis* on the one hand respire the excess oxygen produced by the cyanobacterial farmer, while triggering and improving biofilm formation of laborer and farmer organism. Similar findings have already been reported for a different reaction system.<sup>5</sup>

For the cultivation of the proposed consortium a reactor system had to be developed that meets the consortiums and the reaction needs. Most important is the gas tightness, as our product H<sub>2</sub> is very volatile and escapes easily through many materials and connectors and is only produced under anaerobic conditions. On the other hand, a transfer to outdoor conditions of the reactor system is important for the

---

<sup>1</sup> Flemming, H.-C. & Wuertz, S. Bacteria and archaea on Earth and their abundance in biofilms. *Nat Rev Microbiol* (2019)17: 247

<sup>2</sup> The biofilm matrix: multitasking in a shared space. *Nat Rev Microbiol* (2023) 21:70

<sup>3</sup> Philipp, L.-A., Bühler, K., Ulber, R. & Gescher, J. Beneficial applications of biofilms. *Nat Rev Microbiol* *accepted*

<sup>4</sup> Brown, B., Wilkins, M., Saha, R. *Rhodopseudomonas palustris*: A biotechnology chassis. *Biotechnol Adv* (2022) 60: 108001

<sup>5</sup> Hoschek, A., et al. Mixed-species biofilms for high-cell-density application of *Synechocystis* sp. PCC 6803 in capillary reactors for continuous cyclohexane oxidation to cyclohexanol. *Bioresour Technol* (2019) 282:171

use of sunlight as energy source for a cost-efficient production system. In the following the capillary biofilm reactor (CBR) is described and a standard operation procedure is reported.

## 2. Methodology

### 2.1 Strains used in this study

Species	Strain	References
<i>Rhodopseudomonas palustris</i>	DSM127	Malik et al. (1987) <sup>6</sup>
<i>Pseudomonas taiwanensis</i>	VLB120_eGFP	Heuschkel, Hoschek et al. (2019) <sup>7</sup>

### 2.2 Growth media used

The following media were used for batch cultivation of *R. palustris* (R8AH-H2), batch cultivation of *P. taiwanensis* (M9-Medium) and biofilm cultivation (M27-H2). All media are prepared under sterile conditions.

Table 2: R8AH-H<sub>2</sub> Medium for cultivation of *R. palustris*

Component	Formula	M [g/mol]	Conc. [g/L]	Volume [mL/L]
Sodium acetate	C <sub>2</sub> H <sub>3</sub> NaO <sub>2</sub>	82.03	1	
Yeast extract	C <sub>19</sub> H <sub>14</sub> O <sub>2</sub>	319.19	1	
Magnesium sulfate	MgSO <sub>4</sub> x 7 H <sub>2</sub> O	246.48	0.2	
Calcium chloride	CaCl <sub>2</sub> x 2 H <sub>2</sub> O	147.01	0.07	
Di-sodium EDTA dihydrate	C <sub>10</sub> H <sub>14</sub> N <sub>2</sub> Na <sub>2</sub> O <sub>8</sub> x 2H <sub>2</sub> O	372.21	0.023	
Potassium Dihydrogen Phosphate	KH <sub>2</sub> PO <sub>4</sub>	136.09	0.6	
Di-Potassium hydrogen phosphate	K <sub>2</sub> HPO <sub>4</sub>	174.18	0.9	
Fe(III) citrate solution	C <sub>6</sub> H <sub>5</sub> FeO <sub>7</sub>	244.94	0.01	
R8AH Trace element solution				1.0
R8AH Viatmin solution				7.5

#### R8AH Trace element solution

Iron(III)citrate	C <sub>6</sub> H <sub>5</sub> FeO <sub>7</sub>	244.94	3
Manganese sulfate monohydrate	MnSO <sub>4</sub> x H <sub>2</sub> O	169.02	0.02
Boric acid	H <sub>3</sub> BO <sub>3</sub>	61.83	0.01
Copper sulfate pentahydrate	CuSO <sub>4</sub> x 5 H <sub>2</sub> O	249.68	0.01

<sup>6</sup> Malik, K.A., et al. A survey of hydrogen production, nitrogen fixation and hydrogen metabolism in Rhodospirillaceae. In: Neijssel, O.M. (eds). Proc. 4th European Congress on Biotechnology. (1987) 3: 558; Elsevier, Amsterdam

<sup>7</sup> Heuschkel, I. et al. Data on mixed trophies biofilm for continuous cyclohexane oxidation to cyclohexanol using *Synechocystis* sp. PCC 6803. Data Brief (2019) 25: art. 104059

Ammonium molybdate tetrahydrate	$(\text{NH}_4)_6\text{Mo}_7\text{O}_{24} \times 4 \text{H}_2\text{O}$	1235.86	0.02
Zinc sulfate heptahydrate	$\text{ZnSO}_4 \times 7 \text{H}_2\text{O}$	287.60	0.018
EDTA	$\text{C}_{10}\text{H}_{16}\text{N}_2\text{O}_8$	292.24	0.5
Calcium chloride	$\text{CaCl}_2 \times 2 \text{H}_2\text{O}$	147.01	0.2
<b>R8AH Vitamin solution</b>			
Nicotinic acid (B3)	$\text{C}_6\text{H}_5\text{NO}_2$	123.11	0.2
Nicotinamide	$\text{C}_6\text{H}_6\text{N}_2\text{O}$	122.12	0.2
Thiamin-HCl dihydrate (B1)	$\text{C}_{12}\text{H}_{22}\text{Cl}_2\text{N}_4\text{O}_3\text{S}$	373.30	0.4
Biotin (B7)	$\text{C}_{10}\text{H}_{16}\text{N}_2\text{O}_3\text{S}$	244.31	0.008

Table 3: M27-H<sub>2</sub> Medium for the cultivation of dual species biofilms.

Component	Formula	M [g/mol]	Conc. [g/L]	Volume [mL/L]
Yeast extract	$\text{C}_{19}\text{H}_{14}\text{O}_2$	319.193	0.3	
Sodium acetate	$\text{C}_2\text{H}_3\text{NaO}_2$	82.0343	1	
Potassium dihydrogen phosphate	$\text{KH}_2\text{PO}_4$	136.086	0.5	
magnesium sulfate	$\text{MgSO}_4 \times 7 \text{H}_2\text{O}$	246.48	0.4	
Sodium chloride	$\text{NaCl}$	58.443	0.4	
Calcium chloride	$\text{CaCl}_2 \times 2 \text{H}_2\text{O}$	147.01	0.05	
Fe(III) citrate solution	$\text{C}_6\text{H}_5\text{FeO}_7$	244.944	1	
Vitamin B12-solution	$\text{C}_{63}\text{H}_{88}\text{CoN}_{14}\text{O}_{14}\text{P}$	1355.38	0.1	
L-Cysteiniumchloride	$\text{C}_3\text{H}_{10}\text{ClNO}_3\text{S}$	175.64	300	
Trace element solution SL-6				1
<b>Trace element solution SL-6</b>				
Zinc sulfate heptahydrate	$\text{ZnSO}_4 \times 7 \text{H}_2\text{O}$	287.60	0.1	
Mangan(II)chloride tetrahydrate	$\text{MnCl}_2 \times 4 \text{H}_2\text{O}$	197.91	0.03	
Borate	$\text{H}_3\text{BO}_3$	61.83	0.3	
Cobalt chloride hexahydrate	$\text{CoCl}_2 \times 6 \text{H}_2\text{O}$	237.93	0.2	
Copper chloride dihydrate	$\text{CuCl}_2 \times 2 \text{H}_2\text{O}$	170.48	0.01	
Nickel chloride hexahydrate	$\text{NiCl}_2 \times 6 \text{H}_2\text{O}$	237.69	0.02	
Sodium molybdate	$\text{Na}_2\text{MoO}_4 \times 2 \text{H}_2\text{O}$	241.95	0.03	

For the cultivation in the optimized CBR system, the medium is split into carbon source and salt medium. Both are two times concentrated.

Table 4: M9-Medium for the cultivation of *P. taiwanensis*.

Component	Formula	M [g/mol]	Conc. [g/L]	Volume [mL/L]
-----------	---------	-----------	-------------	---------------

10x M9 salt solution			100
Magnesium sulfate (1M)	MgSO <sub>4</sub>	120.37	2
US* trace elements solution			1
Thiamine (1 % w/v)	C <sub>12</sub> H <sub>17</sub> N <sub>4</sub> OS <sup>+</sup>	265.36	1
Glucose (50% w/v)	C <sub>6</sub> H <sub>12</sub> O <sub>6</sub>	180.16	4
Water	H <sub>2</sub> O	18.02	982

#### 10x M9 salt solution

Disodium hydrogen phosphate heptahydrate	Na <sub>2</sub> HPO <sub>4</sub> x 7 H <sub>2</sub> O	268.07	128.00
Potassium dihydrogen phosphate	KH <sub>2</sub> PO <sub>4</sub>	136.09	30.00
Sodium chloride	NaCl	58.443	5.00
Ammonium chloride	NH <sub>4</sub> Cl	53.49	10.00

#### US\* trace elements solution

Hydrochloric acid (37 % w/v)	HCl	12.02	83.20
Iron (II) sulfate dihydrate	FeSO <sub>4</sub> x 2 H <sub>2</sub> O	187.94	4.87
Calcium chloride dihydrate	CaCl <sub>2</sub> x 2 H <sub>2</sub> O	147.01	4.10
Mangan chloride tetrahydrate	MnCl <sub>2</sub> x 4 H <sub>2</sub> O	197.91	1.50
Zinc sulfate	ZnSO <sub>4</sub>	161.47	1.05
Borate	H <sub>3</sub> BO <sub>3</sub>	201.22	0.30
Disodium molybdate dihydrate	Na <sub>2</sub> MoO <sub>4</sub> x 2 H <sub>2</sub> O	241.95	0.25
Copper chloride dihydrate	CuCl <sub>2</sub> x 2 H <sub>2</sub> O	134.45	0.15
Disodium EDTA dihydrate	Na <sub>2</sub> EDTA x 2 H <sub>2</sub> O	372.24	0.84

### 3. Standard operation procedure CBR:

SOP for the CBR will be described in 3 consecutive steps, the preparation of the inoculum, the reactor assembly and the reactor run.


#### 1. Preparation of the inoculum:

- a. Cultivation of *R. palustris* preculture:
  1. For preculture 1 mL of *R. palustris* cryo culture is inoculated into 50 mL R8AH-H2 Medium in 120 mL serum flask filled with N<sub>2</sub> gas.
  2. Cultivation was performed at 30 °C, 50 μE m<sup>-2</sup> s<sup>-1</sup> for 6 days.
  3. From these precultures inoculation cultures in 10 mL M27-H2 medium with an OD<sub>650</sub>=2 were prepared 12 h before inoculation and cultivated like the preculture.
- b. Cultivation of *P. taiwanensis* preculture:
  1. For the first preculture 5 mL LB medium are inoculated with an inoculation loop of a *P. taiwanensis* preculture and cultivated at 30 °C, 200 rpm (2.5 cm amplitude) over night.

- 10 mL M9 medium are inoculated with 200  $\mu$ L of the first preculture and cultivated at 30  $^{\circ}$ C, 200 rpm (2.5 cm amplitude) for 24 h.
- 10 mL M9 medium are again inoculated with 200  $\mu$ L of step 1.b.2 and cultivated for 6-8 h at 30  $^{\circ}$ C, 200 rpm (2.5 cm amplitude)
- 10 mL M27-H2 medium are inoculated with an  $OD_{450}=2$  12 h before inoculation and cultivated like described for the precultures.

## 2. Reactor assembly:

For running two CBR's the use of two 1L split medium bottles and one 2L waste bottle is recommended. Schematic overview figure 1.



**Figure 1:** Simplified schematic overview of the CBR system for  $H_2$  production using *R. palustris* and *P. taiwanensis*. The CBR cultivation unit is located in a gas tight, transparent glas casket which allows anaerobic cultivation as it is constantly flushed with either  $N_2$  or Ar depending on the cultivation conditions. The gaseous product is collected in the gas bags connected to either sample port 1 or 2.

- The tubings for the reactors are assembled according to Heuschkel et al. 2019<sup>8</sup> with the addition of tubing for the additional medium bottle (2x medium lines, Tygon, 20cm length, Inner diameter: 2mm) and a second t-connector in front of the t-connector that connects the gas peristaltic pump tube (PVC, inner diameter 2.01 mm, Watson Marlow) with the inlet of the CBR that connects the two split medium flows before the reactor inlet.
- The assembly of the reactor unit is performed before autoclaving. The reactor consists of an outer glass tube with two GL45 connections on the ends.
- A GL45 3 port cap is attached to the outlet end of the reactor unit with two GL14 hose connectors with PTFE fitting both equipped with a PTFE tube (PTFE, 5 cm length, Inner diameter: 2mm). One connects to the silicon capillary that is located within the outer glass tube, the other is used as outlet for the glass mantle stream of gas. One GL14 screw cap with a hole is supplied with a silicon membrane coated with PTFE for gas sampling of the mantle space.
- A GL45 2 port cap is attached to the inlet of the glass tube supplied with two GL45 hose connectors with PTFE fitting both equipped with a PTFE tube (PTFE, 5 cm length, Inner diameter: 2mm) as used on the outlet. One is used to connect to the silicon capillary (silicon, length: 22 cm, inner

<sup>8</sup> Heuschkel, I. et al. Data on mixed trophies biofilm for continuous cyclohexane oxidation to cyclohexanol using *Synechocystis* sp. PCC 6803. Data Brief (2019) 25: art. 104059

diameter: 3mm) that is located within the outer glass tube, the other is used as inlet for the glass mantle stream of gas. This PTFE tube is on the outside attached to a peristaltic pump tube (PVC, inner diameter 2.01 mm, Watson Marlow) as used for the medium pumps and a sterile filter. This line can be used to pump any gas or simply air through the mantle of the reactor.

5. After autoclaving, the reactor is assembled under sterile conditions.

### 3. Preparation for operation:

1. The medium bottles are attached to the CBR setup. The pumps are turned on and the entire system is filled with gas and medium. The gas tubes except for mantle flow (constantly on at  $250 \mu\text{L min}^{-1}$ ) are tightly closed using metal clamps and only medium is flushed through the capillary for 1-2 h at a flow rate of  $52 \mu\text{L min}^{-1}$  prior to inoculation.
2. The two inoculation cultures are combined in a luer lock syringe (2 mL of each preculture). A cannula is attached to the luer lock syringe. The safety cap is kept on until inoculation.
3. For inoculation the medium flow is stopped and the 3-way valve in front of the inoculation port is closed. Inoculation takes place next to a flame to ensure sterile conditions. The inoculation port is sterilized with 70 % ethanol solution. The cannulas safety cap is removed and 2 mL of inoculum is injected into each capillary.
4. The CBR's are covered with aluminum foil or any other foil that prohibits light transmission to the capillary and left settling for 24 h.
5. After 24 h the foil is removed and medium flow is started with a final combined flow of  $52 \mu\text{L min}^{-1}$ . If run under indoor conditions the overhead LED's are turned on with  $50 \mu\text{mol m}^{-2} \text{s}^{-1}$ . The gas tubes are kept closed and the reactor is run in continuous flow mode for 72 h.
6. After these 72 h in constant flow, the metal clamps are removed from the gas tubes and the gas flow is started with a flow of  $52 \mu\text{L min}^{-1}$ .

### Sampling options:

1. The gas phase of the reactor mantle and the gas phase of the capillary can be sampled through the membrane in the reactor mantle and the bubble trap at the capillary outlet and either quantified via gas chromatography or membrane inlet mass spectrometry (MIMS). MIMS can additionally be installed in the liquid phase of the reactor to quantify dissolved gasses.
2. The outflow of the reactor can be sampled by the introduction of an additional 3-way valve after the capillary. The cell count and populations can be directly analyzed via flow cytometry. After centrifugation, the cell-free supernatant can be used for quantification of extracellular metabolites, components and salts.
3. During the run, the biofilm thickness can be monitored via optical coherence tomography.
4. After the reactor run, endpoint measurements including omics and EPS analysis of the biofilm can be performed.